

# HOW TO USE IRRIMED:

The application is rather easy because it is similar compared to conventional cleaning/rinsing methods.

Fill your irrimed with nicely tempered tap water, which is then pumped via a transfer tube through a cone into the stoma.

This stimulates a stretch, which in turn activates the peristaltic and therefore also a bowel emptying.

Consequently, the colon has been emptied, allowing you to achieve stool-free periods of 24 to 48 hours and more.

You can find a detailed description of the process in our instruction manual, see:

[www.schaefers-medizintechnik.com](http://www.schaefers-medizintechnik.com)



## OUR SERVICE PLEDGE

- **Competent counselling**
- **Repairs in less than a week**
- **Free replacement units**

### Order details:

irrimed rectal  
Order number: 4000

Irrimed stoma and all needed accessories are available at specialist dealers or for self-paying customers directly with us. Irrimed stoma can be prescribed by prescription.

irrimed classic

A PRODUCT BY:



**Schaefers**  
Medizintechnik GmbH

Im Wiesengrund 13  
46342 Velen  
Telefon: +49(0)2863 95024

service@schaefers-medizintechnik.com  
www.schaefers-medizintechnik.com

As of 2024-02-06

irrimed classic

# STOMA WITHOUT A BAG... THE SOLUTION: IRRIMED



# CLEAN, SAFE, FEEL-GOOD

## IRRIGATION WITH THE IRRIMED

### WHAT IS IRRIGATION?

*Irrigation is the best and most efficient method to regain controlled bowel emptying: It is possible to achieve stool-free periods of 24 to 48 hours or even more.*

*Irrimed prevents typical problems of bag supply, such as*

*... noise nuisance*

*... unpleasant odour*

*... skin irritations due to excreta*

*... disposal problems of the bag etc.*

*the stoma can be covered discreetly with a stoma cap.*

*Irrimed offers more security, comfort and allows you to enjoy your private and social life.*

*Additionally, it has a positive impact on the environment by using fewer resources.*



### INDEPENDENT

Thanks to its modern and powerful Li-Ion battery, irrimed stoma can be used anywhere, anytime, autonomous of the mains supply

### COMFORTABLE

You can easily check the water temperature from the large LCD display.

### COMPACT

While travelling simply remove the handle and safe space.

### MOBILE

The ergonomically optimal water container and separate handle allow you to carry the full container anywhere you like, with ease.

### SECURE

You can fully recharge the battery in max. 8 hours, with the quick charger, which can be used world-wide.